


Preprocessing and extraction workflows for the MagMAX Cell-Free DNA Isolation Kit

Recommended protocols for processing samples prior to nucleic acid extraction via the Applied Biosystems™ MagMAX™ Cell-Free DNA Isolation Kit



* 2 mL of sample input = 35-minute run time on a Thermo Scientific™ KingFisher™ magnetic particle processor.
 4–5 mL of sample input = 40–50-minute run time on a KingFisher magnetic particle processor.

For Research Use Only. Not for use in diagnostic procedures. © 2024 Thermo Fisher Scientific Inc. All rights reserved.
 All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Streck is a trademark of Streck Laboratories Inc.

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk
 Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it
 Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl
 Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es
 Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2024 Thermo Fisher Scientific Inc. All rights reserved.
 Trademarks used are owned as indicated at fishersci.com/trademarks.