

**Single-use technology**

Standard single-use bottle assembly systems

Fluid transfer systems can easily consist of multiple components that all need to be sourced, purchased, inventoried, and assembled. With the Thermo Scientific™ Standard Single-Use Bottle Assembly Systems, we can relieve you of the burden of managing individual components and validating the systems in-house. Realize labor savings and mitigate risk by leaving the assembly and processing to us. These easy-to-order configured solutions include rigid support containers with preconfigured cap sets, and are delivered gamma-irradiated for use right out of the package.

Applications

- Media and buffer transfer
- Product sampling
- Harvest collection and bulk storage
- Filtration and purification process connection

Key benefits

- Designed for a broad range of applications, including bioproduction, life science, and general lab use
- Helps eliminate setup and post-use cleaning steps required for reusable containers
- Helps reduce cross-contamination risks
- Allows process-specific flexibility with multiple bottle sizes

Assembly details

- Bottles are available in polycarbonate (PC) and polyethylene terephthalate glycol (PETG) contact surfaces
- Each bottle is topped with a cap assembly including tubing, plug, and filter
- Most assemblies contain two ports with C-Flex™ tubing; the 2 L size includes a third port, also with C-Flex tubing
- One line on each assembly is dedicated to venting
- The second and third lines are for liquid transfer; one plugged line has a dip tube and includes plug fittings, to be used with your preferred welding or sealing equipment
- Fluid path of the finished product has been validated according to AAMI TIR33 with a 10^{-6} sterility assurance level (SAL)

Ordering information

Size	Material	Closure	Tubing connection	Tubing length, in. (cm)	Fluid path ID x OD, in. (mm)	Fisher Scientific Cat. No.	Cat. No. (case qty.)
125 mL	PC	38–430 mm	Port 1: C-Flex tubing with plug Port 2: Vent filter	Port 1: 24 (60.96) Port 2: 2 (5.08)	Port 1: 1/8 x 1/4 (3.2 x 6.4) Port 2: 1/8 x 1/4 (3.2 x 6.4)	17877371	SB00001-I (25/case)
250 mL		38–430 mm	Port 1: C-Flex tubing with plug Port 2: Vent filter	Port 1: 24 (60.96) Port 2: 2 (5.08)	Port 1: 1/8 x 1/4 (3.2 x 6.4) Port 2: 1/8 x 1/4 (3.2 x 6.4)	17887371	SB00002-I (10/case)
500 mL		38–430 mm	Port 1: C-Flex tubing with plug Port 2: Vent filter	Port 1: 24 (60.96) Port 2: 2 (5.08)	Port 1: 1/8 x 1/4 (3.2 x 6.4) Port 2: 1/8 x 1/4 (3.2 x 6.4)	17897371	SB00003-I (12/case)
1 L		38–430 mm	Port 1: C-Flex tubing with plug Port 2: Vent filter	Port 1: 24 (60.96) Port 2: 2 (5.08)	Port 1: 1/8 x 1/4 (3.2 x 6.4) Port 2: 1/8 x 1/4 (3.2 x 6.4)	17807381	SB00004-I (12/case)
2 L		53 B	Port 1: C-Flex tubing with plug Port 2: Vent filter Port 3: C-Flex tubing with plug	Port 1: 12 (30.48) Port 2: 3 (7.62) Port 3: 12 (30.48)	Port 1: 1/4 x 3/8 (6.4 x 9.6) Port 2: 1/4 x 3/8 (6.4 x 9.6) Port 3: 1/4 x 3/8 (6.4 x 9.6)	17817381	SB00005-I (8/case)
125 mL	PETG	38–430 mm	Port 1: C-Flex tubing with plug Port 2: Vent filter	Port 1: 24 (60.96) Port 2: 2 (5.08)	Port 1: 1/8 x 1/4 (3.2 x 6.4) Port 2: 1/8 x 1/4 (3.2 x 6.4)	17827381	SB00006-I (25/case)
250 mL		38–430 mm	Port 1: C-Flex tubing with plug Port 2: Vent filter	Port 1: 24 (60.96) Port 2: 2 (5.08)	Port 1: 1/8 x 1/4 (3.2 x 6.4) Port 2: 1/8 x 1/4 (3.2 x 6.4)	17837381	SB00007-I (20/case)
500 mL		38–430 mm	Port 1: C-Flex tubing with plug Port 2: Vent filter	Port 1: 24 (60.96) Port 2: 2 (5.08)	Port 1: 1/8 x 1/4 (3.2 x 6.4) Port 2: 1/8 x 1/4 (3.2 x 6.4)	17847381	SB00008-I (10/case)
1 L		38–430 mm	Port 1: C-Flex tubing with plug Port 2: Vent filter	Port 1: 24 (60.96) Port 2: 2 (5.08)	Port 1: 1/8 x 1/4 (3.2 x 6.4) Port 2: 1/8 x 1/4 (3.2 x 6.4)	17857381	SB00009-I (20/case)
2 L		53 B	Port 1: C-Flex tubing with plug Port 2: Vent filter Port 3: C-Flex tubing with plug	Port 1: 12 (30.48) Port 2: 3 (7.62) Port 3: 12 (30.48)	Port 1: 1/4 x 3/8 (6.4 x 9.6) Port 2: 1/4 x 3/8 (6.4 x 9.6) Port 3: 1/4 x 3/8 (6.4 x 9.6)	17867381	SB00010-I (8/case)

Learn more at eu.fishersci.com

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk
Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it
Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl
Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es
Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

