



Bioprocessing

3D Productainer BioProcess Containers

A proven solution leading the industry in quality and reliability

3D Thermo Scientific™ Productainer™ BioProcess Containers (BPCs) are used to confidently store, transport, and ship media, buffers, process liquids, and biological products. With design flexibility throughout the comprehensive portfolio of these BPCs, you can trust our high-quality films to provide robust and reliable liquid containment solutions as you scale from 50 L to 1,000 L.

Quality:







- Proven, high-quality film featuring leading single-use Thermo Scientific™ Aegis™ 5-14 film and Thermo Scientific™ CX5-14 film
- Films optimized to enable excellent and reproducible cell growth
- High level of assurance with automated manufacturing and lot-based endotoxin and particulate testing
- Minimized risk of leaks and simplified setup with the Thermo Scientific™ BioTitan™ Retention Device, which provides a 360° seal

Reliability:

- Versatile range of standard sizes that are designed to fit a full range of support containers
- Strength, flexibility, and reliability for media, buffer, harvest, drug substance, and drug product processes
- High-quality products that are ready to ship
- Leading single-use manufacturing footprint for redundancy and supply assurance

Find it at eu.fishersci.com

Ordering information

	Description	Cat. No	Fisher Scientific Cat. No
	Productainer BioProcess Container (BPC), 200 L with 2 Top and 1 Bottom Ports, 3D Single-Use, CX5-14 film	SH30650.03	15556829
	Productainer BioProcess Container (BPC), 100 L with 2 Top and 1 Bottom Ports, 3D Single-Use, CX5-14 film	SH30650.02	15546829
	Productainer BioProcess Container (BPC), 200 L with 2 Top and 1 Bottom Ports (with tri-clamps), 3D Single-Use, CX5-14 film	SH30672.03	15506849
	Productainer BioProcess Container (BPC), 200 L with 4 Top Ports and 1 Bottom Dispense, 3D Single-Use, CX5-14 film	SH30652.03	15536839
	Productainer BioProcess Container (BPC), 200 L with 2 Top Ports, 3D Single-Use, CX5-14 film, single layer	SH30649.03	15526829
	Productainer BioProcess Container (BPC), 200 L with 4 Top Ports, 3D Single-Use, CX5-14 film, single layer	SH30651.03	15586829

Learn more at eu.fishersci.com

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at **Belgium:** fishersci.be **Denmark:** fishersci.dk
Germany: fishersci.de **Ireland:** fishersci.ie **Italy:** fishersci.it
Finland: fishersci.fi **France:** fishersci.fr **Netherlands:** fishersci.nl
Norway: fishersci.no **Portugal:** fishersci.pt **Spain:** fishersci.es
Sweden: fishersci.se **Switzerland:** fishersci.ch **UK:** fishersci.co.uk

 **fisher scientific**
part of Thermo Fisher Scientific